



Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support



100% made in Germany

TC Axial Handle Endoscopic Needle Holder


Main +1 (800) 303-3050 Fax +1 (949) 600-8874 9620 Research Drive, Irvine, CA 92618 USA
info@sbhsurgical.com www.sbhsurgical.com

Copyright © SBH ®. 2020 - All rights reserved.

TC NEEDLE HOLDER WITH AXIAL RATCHET HANDLE

Lightweight and reliable Needle Holders with two different easy rotatable handles for a maximum of flexibility and comfort. These Needle Holders are designed for a precise handling to grasp and hold needles more securely. The flush port ensures an easy cleaning.

Ø 5mm
330mm working length
Luer-Lock connection



Ergonomic, pistol grip handle design for improved comfort and effortless working

Straight handle design: ratchet can be opened either with thumb or forefinger.



300-300-050

300-300-000

TC inserts, cross serrated
STRAIGHT



300-300-150

300-300-100

TC inserts, cross serrated
LEFT CURVED



300-300-250

300-300-200

TC inserts, cross serrated
RIGHT CURVED



300-300-550

300-300-500

Needle Positioner, serrated,
convex-concave jaws, diamond dusted
STRAIGHT