



Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support



SBH 100-120 Ring Applicator for Tubal Sterilization, for single puncture – for use with Operating Laparoscopes with \varnothing 6 mm instrument channels (SBH 099-006 and SBH 099-106). A knurled screw is fixing the Ring Applicator securely to the telescope.

SBH 100-122 Ring Applicator for Tubal Sterilization, for second puncture – for use with Trocar sleeves \varnothing 7.0 mm

SBH 100-124 Plastic Rings, package of 50 pieces, autoclavable



SBH 100-126 Plastic Rings, sterile package of 2 pieces

All Ring Applicators are supplied with a Quick Loading Adaptor to ensure a fast and easy loading of the applicator.

100% made in Germany

Fallopian Ring Applicator

Main +1 (800) 303-3050 Fax +1 (949) 600-8874 9620 Research Drive, Irvine - CA 92618 USA info@sbhsurgical.com www.sbhsurgical.com

Copyright © SBH ®. 2020- All rights reserved.