



Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support

- Ergonomic
- Lightweight Aluminum and Heavy Axial Handle



100% made in Germany

Laparoscopic Needle Holders


Main: +1 (800) 303-3050 Fax: +1 (949) 600-8874 9620 Research Drive, Irvine, CA 92618 USA
info@sbhsurgical.com www.sbhsurgical.com

Copyright © SBH ®. 2021 - All rights reserved.




Needleholder "non detachable"


Needleholder with axial handle for the area of the Thoracic and Cardio Vascular Surgery



On request - Axial handle also without ratchet as well as in different colors available



Luer Lock flushing port for optimal cleaning




ratchet for always optimal positioning, with overload protection




round handle for more comfortable grip easy handling



Description	∅	length	complete instrument
tip straight, TC	5 mm	300 mm	SBH 30-50440.30
	3 mm	300 mm	SBH 30-30440.30
tip curved left, TC	5 mm	300 mm	SBH 30-50441.30
	3 mm	300 mm	SBH 30-30441.30
tip curved right, TC	5 mm	300 mm	SBH 30-50442.30
	3 mm	300 mm	SBH 30-30442.30



Needleholder "non detachable"



Needleholder with axial handle grey - on request



Luer Lock flushing port for optimal cleaning




toothed ratchet for always optimal positioning, with overload protection





aluminium anodized axial handle

- Alu-handle looks heavy model but light weight
- Easy handling more comfort with overload protection







Needleholder "non detachable"



Needleholder, Orca with axial handle



aluminium anodized handle, straight, for more comfortable grip allow easy handling



aluminium anodized handle, curved, for more comfortable grip allow easy handling



Luer Lock flushing port for optimal cleaning



toothed ratchet for always optimal positioning, with overload protection













tip with 1x2 teeth additional security

Needleholder "non detachable"



Needleholder with heavy axial handle



Description	Ø	length	complete instrument	complete instrument	complete instrument
 tip straight, TC	5 mm	330 mm	SBH 30-50400.33	SBH 30-50410.33	SBH 30-50143.33
		450 mm	SBH 30-50400.45	SBH 30-50410.45	SBH 30-50143.45
 tip straight, TC	3 mm	220 mm	SBH 30-30400.22	SBH 30-30410.22	SBH 30-30143.22
		330 mm	SBH 30-30400.33	SBH 30-30410.33	SBH 30-30143.33
 tip curved left, TC	5 mm	330 mm	SBH 30-50401.33	SBH 30-50411.33	SBH 30-50142.33
		450 mm	SBH 30-50401.45	SBH 30-50411.45	SBH 30-50142.45
 tip curved left, TC	3 mm	220 mm	SBH 30-30401.22	SBH 30-30411.22	SBH 30-30142.22
		330 mm	SBH 30-30401.33	SBH 30-30411.33	SBH 30-13042.33
 tip curved right, TC	5 mm	330 mm	SBH 30-50402.33	SBH 30-50412.33	SBH 30-50188.33
		450 mm	SBH 30-50402.45	SBH 30-50412.45	SBH 30-50188.45
 tip curved right, TC	3 mm	220 mm	SBH 30-30402.22	SBH 30-30412.22	SBH 30-30188.22
		330 mm	SBH 30-30402.33	SBH 30-30412.33	SBH 30-30188.33
 Flamingo tip straight	5 mm	330 mm	SBH 30-50403.33	SBH 30-50413.33	SBH 30-50189.33
		450 mm	SBH 30-50403.45	SBH 30-50413.45	SBH 30-50189.45
 tip straight 1x2 teeth, TC	5 mm	330 mm	SBH 30-50404.33	SBH 30-50414.33	SBH 30-50241.33
		450 mm	SBH 30-50404.45	SBH 30-50414.45	SBH 30-50241.45
 tip straight, serrated	5 mm	330 mm	SBH 30-50405.33	SBH 30-50415.33	SBH 30-50187.33
		450 mm	SBH 30-50405.45	SBH 30-50415.45	SBH 30-50187.45




Needleholder 3-part detachable






SBH 30-5300-033
Metal shaft Ø 5 mm, 330 mm length

SBH 30-5300-045
Metal shaft Ø 5 mm, 450 mm length




Description	Ø	length
tip straight, TC	5 mm	330 mm
		450 mm
tip curved left, TC	5 mm	330 mm
		450 mm
tip curved right, TC	5 mm	330 mm
		450 mm
tip straight „Flamingo“	5 mm	330 mm
		450 mm
tip straight 1x2 teeth, TC	5 mm	330 mm
		450 mm
tip straight, serrated	5 mm	330 mm
		450 mm


Needleholder 3-part detachable



SBH 30-4001Z
Axial handle, Orca, straight



SBH 30-4002Z
Axial handle, Orca, curved




SBH 30-4000Z
Heavy axial handle













insert only	complete instrument	complete instrument	complete instrument
SBH 30-53420.33	SBH 30-53400.33	SBH 30-53410.33	SBH 30-53143.33
SBH 30-53420.45	SBH 30-53400.45	SBH 30-53410.45	SBH 30-53143.45
SBH 30-53421.33	SBH 30-53401.33	SBH 30-53411.33	SBH 30-53142.33
SBH 30-53421.45	SBH 30-53401.45	SBH 30-53411.45	SBH 30-53142.45
SBH 30-53422.33	SBH 30-53402.33	SBH 30-53412.33	SBH 30-53188.33
SBH 30-53422.45	SBH 30-53402.45	SBH 30-53412.45	SBH 30-53188.45
SBH 30-53423.33	SBH 30-53403.33	SBH 30-53413.33	SBH 30-53189.33
SBH 30-53423.45	SBH 30-53403.45	SBH 30-53413.45	SBH 30-53189.45
SBH 30-53424.33	SBH 30-53404.33	SBH 30-53414.33	SBH 30-53241.33
SBH 30-53424.45	SBH 30-53404.45	SBH 30-53414.45	SBH 30-53241.45
SBH 30-53425.33	SBH 30-53405.33	SBH 30-53415.33	SBH 30-53187.33
SBH 30-53425.45	SBH 30-53405.45	SBH 30-53415.45	SBH 30-53187.45



Needleholder 2-part detachable



Description	Ø	length	shaft / insert	compl. Instrument	compl. Instrument	compl. Instrument
 tip straight, TC	5 mm	330 mm	SBH 30-52420.33	SBH 30-52400.33	SBH 30-52410.33	SBH 30-52143.33
		450 mm	SBH 30-52420.45	SBH 30-52400.45	SBH 30-52410.45	SBH 30-52143.45
 tip straight, TC	3 mm	220 mm	SBH 30-32420.22	SBH 30-32400.22	SBH 30-32410.22	SBH 30-32143.22
		330 mm	SBH 30-32420.33	SBH 30-32400.33	SBH 30-32410.33	SBH 30-32143.33
 tip curved left, TC	5 mm	330 mm	SBH 30-52421.33	SBH 30-52401.33	SBH 30-52411.33	SBH 30-52142.33
		450 mm	SBH 30-52421.45	SBH 30-52401.45	SBH 30-52411.45	SBH 30-52142.45
 tip curved left, TC	3 mm	220 mm	SBH 30-32421.22	SBH 30-32401.22	SBH 30-32411.22	SBH 30-32142.22
		330 mm	SBH 30-32421.33	SBH 30-32401.33	SBH 30-32411.33	SBH 30-32142.33
 tip curved right, TC	5 mm	330 mm	SBH 30-52422.33	SBH 30-52402.33	SBH 30-52412.33	SBH 30-52188.33
		450 mm	SBH 30-52422.45	SBH 30-52402.45	SBH 30-52412.45	SBH 30-52188.45
 tip curved right, TC	3 mm	220 mm	SBH 30-32422.22	SBH 30-32402.22	SBH 30-32412.22	SBH 30-32188.22
		330 mm	SBH 30-32422.33	SBH 30-32402.33	SBH 30-32412.33	SBH 30-32188.33
 tip straight, Flamingo	5 mm	330 mm	SBH 30-52423.33	SBH 30-52403.33	SBH 30-52413.33	SBH 30-52189.33
		450 mm	SBH 30-52423.45	SBH 30-52403.45	SBH 30-52413.45	SBH 30-52189.45
 tip straight 1x2 teeth, TC	5 mm	330 mm	SBH 30-52424.33	SBH 30-52404.33	SBH 30-52414.33	SBH 30-52241.33
		450 mm	SBH 30-52424.45	SBH 30-52404.45	SBH 30-52414.45	SBH 30-52241.45
 tip straight, serrated	5 mm	330 mm	SBH 30-52425.33	SBH 30-52405.33	SBH 30-52415.33	SBH 30-52187.33
		450 mm	SBH 30-52425.45	SBH 30-52405.45	SBH 30-52415.45	SBH 30-52187.45





Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support

Main +1 (800) 303-3050 Fax +1 (949) 600-8874 9620 Research Drive, Irvine, CA 92618 USA info@sbhsurgical.com www.sbhsurgical.com