



Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support



Instruments for

*Neonatal & Pediatric Surgery
Specialized Newborn Baby's Surgery*

100% made in Germany

Neonatal-Pediatric



Usage and fields of application: PAEDIATRIC SURGERY: NEONATES

Special features:

In this brochure you find a concentration of small/delicate surgical instruments in order to maximise the best outcomes for babies. for quality and efficiency - when it comes to infant and baby care response. Easy to clean while protecting your surgical instruments.

The problem it solves:

- Advances in technology and health care expertise have led to increasing survival rates of very premature babies.
- Clamps available in the two most prevalent atraumatic jaw formations: Cooley and Debakey, providing a secure hold of the vessel without damaging the vessel's delicate intimal lining.



CASTANEDA jaw/tip

- Two rows of microscopic teeth to avoid injury to the vessel



COOLEY jaw/tip

- Two rows of microscopic teeth
- Tissue expands into groove between rows which supports firm hold



DEBAKEY jaw/tip

- Has 1x2 rows of microscopic teeth
- Tissue is held by teeth for a firm grasp



Designed specifically to meet
***the needs of Neonates - when it comes
to infant and baby surgical techniques.***



NEONATAL & PEDIATRIC SURGERY


INDEX

<i>Instrumentation Supplies</i>	4
Baby Mixter Clamp, Metzenbaum Fino Scissors, Forceps, Micro Bulldog, Weitlaner Retractor, Needle Holder	
Weldon Mini Bulldogs	5
Cooley Neonatal Ultralight Clamps	6
Cooley Selman Blood Vessel Clamps	7
Cooley Pediatric Clamps	8
Debakey Neonatal Ring Handle Bulldog	9
Debakey Titanium Neonatal Clamps	10
Debakey Pediatric Clamps	11
Castaneda Neonatal Clamps	12
Castaneda Pediatric Clamps	13
Micro Thorax Retractors	14-15

To meet your specific needs modifications are available upon request.


Common modifications include:

jaw/tip modifications, ratchet modifications, length modifications, special etching





Neonatal-Pediatric



*Baby Mixer Clamp
curved jaws*

SBH RA 9494

140mm



*Metzenbaum Fino
Scissors, curved jaws*

SBH RA 8000

145mm



*Titanium Forceps
1.7mm tip*

SBH RT 0800

150mm



Micro Bulldog

SBH RC 2068

27.2mm



*Weitlaner 2x3
Retractor, sharp teeth*

SBH RS 1500

110mm



**Ryder Mini
TC-Needle Holder**

SBH RN 1810

150mm



**Derf
TC-Needle Holder**

SBH RN 1798

125mm



**Brown
Needle Holder**

SBH RN 1526

130mm



**Brown
Needle Holder**


SBH RN 1550

130mm

WELDON (GLOVER) MINI BULLDOGS




Neonatal-Pediatric



jaw straight

SBH RC 2069


30mm



*jaw sligthly
curved*

SBH RC 2071


28mm



jaw curved

SBH RC 2073


28mm



jaw straight

SBH RC 2076


33mm



*jaw slightly
curved*

SBH RC 2077

31mm



jaw curved

SBH RC 2078


30mm



jaw cvd. back

SBH RC 2079

25mm




jaw sat.

SBH RC 2086

25mm

Available in Titanium: **SBH RT 2076**


SBH RT 2078



jaw 25°

SBH RC 2074


32mm



jaw 45°

SBH RC 2080


33mm



jaw 60°

SBH RC 2081

25mm



jaw 90°

SBH RC 2082


19mm



Container


SBH RE 1010

ca. 300x200x50mm





COOLEY NEONATAL ULTRALIGHT CLAMPS



100%



jaw straight

jaw straight

jaw slightly curved

jaw curved

s-type

handles straight

handles curved

handles straight

handles straight

handles curved

SBH RC 6300

SBH RC 6302

SBH RC 6303

SBH RC 6306

SBH RC 6307

125mm


125mm

125mm

115mm

125mm

100%



jaw 45°

jaw 45°

jaw 90°

jaw satinsky

jaw double angled

handles straight

handles curved

handles straight

handles curved

handles curved, external carotid clamp

SBH RC 6301

SBH RC 6304

SBH RC 6305

SBH RC 6308

SBH RC 6310

120mm

115mm


100mm

115mm


125mm

To meet your specific needs other models are available upon request.

SELMAN COOLEY BLOOD VESSEL NEONATAL CLAMPS



100%



jaw straight

jaw 50°

jaw fully curved

handles straight

handles straight

handles straight

SBH RC 6150

SBH RC 6151




SBH RC 6170

125mm

120mm

120mm

100%



*jaw 90°
2 scale lines*

*jaw 90°
4 scale lines*

*jaw 90°
6 scale lines*

handles straight

handles straight

handles straight

SBH RC 6355

SBH RC 6356

SBH RC 6357

115 / 20 mm

115 / 30mm

115 / 40mm




Neonatal-Pediatric

COOLEY LEFT & RIGHT ANGLED CLAMPS



*Cooley Kitzmiller
Pediatric Clamp*



jaw left jaw right

handles straight handles straight

SBH RC 6622 SBH RC 6620

130mm 130mm



*Cooley Martinez
left*

*Cooley Martinez
right*

handles straight

handles straight

SBH RC 6612

SBH RC 6610


180mm

180mm

COOLEY PEDIATRIC CLAMPS



100%




*jaw straight
5 scale lines*

handles straight

SBH RC 6350

145mm



*jaw 30°
5 scale lines*

handles straight

SBH RC 6351

140mm




*jaw spoon-shaped
5 scale lines*

handles straight

SBH RC 6353

140mm



*jaw satinsky
6 scale lines*

handles curved

SBH RC 6354

136mm

100%



*jaw
straight*

handles straight

SBH RC 6571


142mm



jaw 15°

SBH RC 6572


141mm



jaw 30°

SBH RC 6573


138mm



jaw 55°

SBH RC 6574


138mm



*jaw
curved L*

SBH RC 6578

135mm



*jaw
curved S*


SBH RC 6579

131mm

DEBAKEY RING HANDLE BULI DELICATE JAW TIP: Ø 1.



100%



jaw straight

jaw curved

handles straight

handles straight

handles curved


SBH RC 2182 SBH RC 2186 SBH RC 2195

134mm

130mm

130mm

100%



jaw 90°

jaw 45°

handles straight


handles curved

SBH RC 2192

SBH RC 2189

105mm

122mm








Neonatal-Pediatric

DEBAKEY NEONATAL ULTRALIGHT CLAMPS TITANIUM DELICATE JAW TIP: Ø 1.2 mm




100%



A
handles
angled

B
handles
curved


C
handles
straight



DEBAKEY PEDIATRIC CLAMP



100%




*satinsky
anastomosis jaw
8 scale lines*

handles curved

SBH RC 2196

145 / 14 mm




*satinsky
anastomosis jaw
9 scale lines*

handles curved

SBH RC 2197

145 / 18 mm




*satinsky
anastomosis jaw
10 scale lines*

handles curved

SBH RC 2198

145 / 23 mm



*strongly curved
anastomosis jaw
9 scale lines*


handles curved

SBH RC 2199

145 / 25 mm



100%




DeBakey Mini Beck

handles straight

SBH RC 3750

165 mm




Derra DeBakey

handles curved

SBH RC 3135

170 mm




Derra DeBakey

handles curved

SBH RC 3136

172mm




Derra DeBakey

handles curved


SBH RC 3137

172mm






CASTANEDA Neonatal CLAMPS




100%



handles curved	SBH RC 6360	SBH RC 6362	SBH RC 6316	SBH RC 6363	SBH RC 6314
	120mm	115mm	117mm	105mm	99mm

Neonatal Castaneda miniature instruments for cardiovascular procedures on the neonate and for peripheral vascular surgery. Very fine and delicate jaws minimize pressure on the arterial wall. Jaws are calibrated at 5mm intervals to determine the size of the anastomosis.

100%




handles curved	SBH RC 6364	SBH RC 6315	SBH RC 6312	SBH RC 6317	SBH RC 6318	SBH RC 6319
	115mm	121mm	120mm	115mm	125mm	120mm

CASTANEDA Pediatric CLAMPS



100%



jaw 36°

jaw 60°

handles straight

handles straight


SBH RC 6412

SBH RC 6433

161mm

146mm

100%



*jaw
satinsky S*

*jaw
satinsky M*

*jaw
satinsky L*

handles curved

handles curved

handles curved

SBH RC 6460


SBH RC 6437

SBH RC 6461

152mm


152mm

152mm






Neonatal-Pediatric






Finochietto Micro Neonatal thorax spreader

Art.No.	overall length	overall height	max. spread	blade depth	blade width
SBH RS 3050	59mm	37mm	32mm	10mm	7mm



Finochietto Micro Neonatal thorax spreader


Art.No.	overall length	overall height	max. spread	blade depth	blade width
SBH RS 3077	109mm	70mm	65mm	26mm	12mm



Titanium Mini Thorax Spreader

Art.No.	overall length	overall height	max. spread
Spreader SBH RS 3160	76.5mm	53mm	39mm

SBH RS 3160SET
consisting of:
SBH RS 3160 - spreader frame
SBH RS 3165 SET - 2 blades
SBH RS 3167 SET - 2 blades fenstered





Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support

Main +1 (800) 303-3050 Fax +1 (949) 600-8874 9620 Research Drive, Irvine - CA 92618 USA info@sbhsurgical.com www.sbhsurgical.com

Copyright © SBH ®. 2020 - All rights reserved.