





Strategic Business Holdings, Inc.


...Quality Products, Superior Customer Service, & Support



Top



Bottom



Middle

100% made in Germany

SL Spine Retractor and Distraction Set for Anterior Cervical Fusion (ACF)



Your address for surgical quality instruments
"100% made in Germany".


SBH stands for highest quality, speed
and efficiency. We will find the technically
and economically optimal solution for you.

Get to know the complete
SBH product range!

For further information please visit us at:
www.sbhsurgical.com
or call us directly at: 1-800-303-3050.

Table of Contents:

1 SL-Set (ACF)	4
2 SL ACF Componets	4-9
3 SL ACF Set Configuration	10






SL Set (ACF)

SBH 110-850-00


SL - Set (ACF), Retractor and Distraction Set with Tray

The new SL Set: Spine Retractor and Distraction Set for Anterior Cervical Fusion (ACF).


- Optimal exposure due to use of longitudinal and transverse double hinged retractors.
- Allows unilateral and bilateral retraction in one set.
- Color coded blades for ease of use together with colored depth gauge.
- All instruments are stored in 3 inserts with special Rilsan coating for the coloured blades to protect instrument surface during transport and sterilization.
- Retractors stainless steel and variety of blades in 5 different designs out of aluminum.
- All components also available separately.



Top



Middle



Bottom



Spinal Distractor

Art. No.	Description	Length Tooth Bar	Opening Width
SBH 110-741-15	Caspar right, short	110mm	100mm
SBH 110-741-20	Caspar right, long	130mm	120mm
SBH 110-741-00	Caspar left, long	110mm	100mm
SBH 110-741-05	Caspar left, long	130mm	120mm



Retractor

Transverse retractor, double hinged, lateral feed

Art. No.		Opening Width
SBH 110-850-10		120mm
SBH 110-850-12		150mm
SBH 110-850-14		200mm




SBH 110-850-20

SL Locking handle for baldes, lateral feed




SBH 110-850-05

Longitudinal Retractor Frame lateral feed





SL ACF Components



Drill Guide with black Radel (PPSU)

Art. No.	Description
SBH 110-744-22	Drill guide left
SBH 110-744-24	Drill guide right



SBH 110-746-23

Caspar screwdriver with black Radel handle (PPSU)



Caspar Twist Drill

Article Number	Description	Total Length	Drill Length	Shaft
SBH 110-748-02	Caspar twist drill, long	148 mm	14 mm	8 mm



Caliper

Article Number	Total Length	Depth	Measuring Range
SBH 120-028-20	240 mm	100 mm	0 - 160 mm




Bone Tamps

Article Number	Length	Diameter
SBH 110-851-03	216 mm	3 mm
SBH 110-851-05	216 mm	3 mm
SBH 110-851-08	216 mm	3 mm




SBH 110-850-30

Bone graft harvesting forceps



SBH 110-850-40

Bone graft holder, total length 216 mm, spike 10 mm




SBH 110-860-60

Color coded depth gauge





SBH 20968-94

Tray for SL-Set (ACF), with 3 inserts and lid (without instruments) dimensions in mm: 540 x 253 x 170





SL ACF Componets



SL - Tissue blade, blunt, 23mm wide

Art. No.	Depth	Color Coded	Set
SBH 110-851-30	30		2
SBH 110-851-35	35		2
SBH 110-851-40	40		2
SBH 110-851-45	45		2
SBH 110-851-50	50		2
SBH 110-851-55	55		2
SBH 110-851-60	60		2
SBH 110-851-65	65		2
SBH 110-851-70	70		2
SBH 110-851-75	75		2

SL - Lateral blade, short teeth, 23mm wide

Art. No.	Depth	Color	Set
SBH 110-852-30	30		1
SBH 110-852-35	35		1
SBH 110-852-40	40		1
SBH 110-852-45	45		1
SBH 110-852-50	50		1
SBH 110-852-55	55		1
SBH 110-852-60	60		1
SBH 110-852-65	65		1

SL - Medial Muscle blade, long teeth, 23mm wide

Art. No.	Depth	Color	Set
SBH 110-853-30	30		1
SBH 110-853-35	35		1
SBH 110-853-40	40		1
SBH 110-853-45	45		1
SBH 110-853-50	50		1
SBH 110-853-55	55		1
SBH 110-853-60	60		1
SBH 110-853-65	65		1




SL - Narrow Tissue blade, blunt, 14 mm wide

Art. No.	Depth	Color	Set
SBH 110-854-35	35		2
SBH 110-854-40	40		2
SBH 110-854-45	45		2
SBH 110-854-50	50		2
SBH 110-854-55	55		2
SBH 110-854-60	60		2



SL - Narrow Muscle blade, short teeth, 14 mm wide

Art. No.	Depth	Color	Set
SBH 110-855-35	35		2
SBH 110-855-40	40		2
SBH 110-855-45	45		2
SBH 110-855-50	50		2
SBH 110-855-55	55		2
SBH 110-855-60	60		2





SSL ACF Componets

SBH 110-850-00

SL - Set (ACF), Retractor and Distraction Set with Tray


Category	Article Number	Quantity
Retractor	SBH 110-850-05	1
	SBH 110-850-10	1
	SBH 110-741-15	1
	SBH 110-741-20	1
	SBH 110-741-00	1
	SBH 110-741-05	1
	SBH 110-850-20	2
	SBH 110-744-22	1
	SBH 110-744-24	1
	SBH 110-746-23	2
	SBH 110-748-02	2
Instruments	SBH 110-028-20	1
	SBH 110-851-03	1
	SBH 110-851-05	1
	SBH 110-850-08	1
	SBH 110-850-30	1
	SBH 110-860-40	1
	SBH 110-860-60	1
Blades	SBH 110-851-30	2
	SBH 110-851-35	2
	SBH 110-851-40	2
	SBH 110-851-45	2
	SBH 110-851-50	2
	SBH 110-851-55	2
	SBH 110-851-60	2
	SBH 110-851-65	2
	SBH 110-851-70	2
	SBH 110-851-75	2
	SBH 110-852-30	1

Category	Article Number	Quantity
Blades	SBH 110-852-35	1
	SBH 110-852-40	1
	SBH 110-852-45	1
	SBH 110-852-50	1
	SBH 110-852-55	1
	SBH 110-852-60	1
	SBH 110-852-65	1
	SBH 110-853-30	1
	SBH 110-853-35	1
	SBH 110-853-40	1
	SBH 110-853-45	1
	SBH 110-853-50	1
	SBH 110-853-55	1
	SBH 110-853-60	1
	SBH 110-853-65	1
	SBH 110-854-35	2
	SBH 110-854-40	2
	SBH 110-854-45	2
	SBH 110-854-50	2
	SBH 110-854-55	2
Tray	SBH 110-855-35	2
	SBH 110-855-40	2
	SBH 110-855-45	2
	SBH 110-855-50	2
	SBH 110-855-55	2
	SBH 110-855-60	2
	SBH 120-968-94	1



Notes:

Handwritten notes area with horizontal lines.





V797D-70161



Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support

Main +1 (800) 303-3050 Fax +1 (949) 600-8874 9620 Research Drive, Irvine, CA 92618 USA info@sbhsurgical.com www.sbhsurgical.com

Copyright © SBH ®. 2021 - All rights reserved.