




Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support



*Interchangeable
with conical
post systems*



*Innovative
Blade-Alignment Groove*

Caspar Cervical Retractor Clear-Line TL



Caspar Cervical Retractor - Clear-Line TL

SBH - 4001 TL Complete



SB-0101	Auto-Lock Transverse Retractor (Standard)	1	SB-2370B	70mm Blunt Blade	2	SB-4655S	55mm ST Multi-Level Blade	2
SB-0104	Auto-Lock Step-down Retractor (Standard)	1	SB-2375B	75mm Blunt Blade	2	SB-4660S	60mm ST Multi-Level Blade	2
SB-0087	Reversible Distractor (LEFT/RIGHT)	1	SB-2330N	30mm Notched Blade	2	SB-4665S	65mm ST Multi-Level Blade	2
SB-0106	Screwdriver	2	SB-2335N	35mm Notched Blade	2	SB-2055	Longitudinal Blade	2
SB-0056	Depth Gage	1	SB-2340N	40mm Notched Blade	2	295-324	Bruns Curette, STR 4-0	1
SB-0082	Bone Tamp	1	SB-2345N	45mm Notched Blade	2	295-333	Bruns Curette, STR 3-0	1
SB-0095	Quick-Release Retractor Handle	2	SB-2350N	50mm Notched Blade	2	295-342	Bruns Curette, STR 2-0	1
SB-1010	Clear-Line Standard Sterilization Container	1	SB-2355N	55mm Notched Blade	2	295-350	Bruns Curette, STR 0	1
SB-0213	Clear-Line - Cloward	2	SB-2360N	60mm Notched Blade	2	295-364	Bruns Curette, ANG 4-0	1
SB-2330B	30mm Blunt Blade	2	SB-2365N	65mm Notched Blade	2	295-373	Bruns Curette, ANG 3-0	1
SB-2335B	35mm Blunt Blade	2	SB-2370N	70mm Notched Blade	2	295-382	Bruns Curette, ANG 2-0	1
SB-2340B	40mm Blunt Blade	2	SB-4630S	30mm ST Multi-Level Blade	2	295-390	Bruns Curette, ANG 0	1
SB-2345B	45mm Blunt Blade	2	SB-4635S	35mm ST Multi-Level Blade	2	295-397	Bruns Curette, REV ANG 4-0	1
SB-2350B	50mm Blunt Blade	2	SB-4640S	40mm ST Multi-Level Blade	2	295-398	Bruns Curette, REV ANG 3-0	1
SB-2355B	55mm Blunt Blade	2	SB-4645S	45mm ST Multi-Level Blade	2	295-399	Bruns Curette, REV ANG 2-0	1
SB-2360B	60mm Blunt Blade	2	SB-4650S	50mm ST Multi-Level Blade	2	295-400	Bruns Curette, REV ANG 0	1
SB-2365B	65mm Blunt Blade	2						



SBH-2004 TL MINI COMPLETE

SB-0101	Auto-Lock Transverse Retractor (Standard)	1
SB-0104	Auto-Lock Step-down Retractor (Standard)	1
SB-0087	Reversible Distractor (LEFT/RIGHT)	1
SB-0106	Screwdriver	2
SB-0056	Depth Gage	1
SB-0082	Bone Tamp	1
SB-0095	Quick-Release Retractor Handle	2
SB-2010	Clear-Line MINI Sterilization Container	1
SB-2330B	30mm Blunt Blade	2
SB-2340B	40mm Blunt Blade	2
SB-2350B	50mm Blunt Blade	2
SB-2360B	60mm Blunt Blade	2
SB-2370B	70mm Blunt Blade	2
SB-2330N	30mm Notched Blade	2
SB-2340N	40mm Notched Blade	2
SB-2350N	50mm Notched Blade	2
SB-2360N	60mm Notched Blade	2
SB-2370N	70mm Notched Blade	2
SB-4630S	30mm ST Multi-Level Blade	2
SB-4635S	35mm ST Multi-Level Blade	2
SB-4640S	40mm ST Multi-Level Blade	2
SB-4645S	45mm ST Multi-Level Blade	2
SB-2055	Longitudinal Blade	2
295-324	Bruns Curette, STR 4-0	1
295-333	Bruns Curette, STR 3-0	1
295-364	Bruns Curette, ANG 4-0	1
295-373	Bruns Curette, ANG 3-0	1
295-397	Bruns Curette, REV ANG 4-0	1
295-398	Bruns Curette REV ANG 3-0	1



SBH-4003 - TL-X

SB-0101	Auto-Lock Transverse Retractor (Standard)	1
SB-0104	Auto-Lock Step-down Retractor (Standard)	1
SB-0056	Depth Gage	1
SB-0095	Quick-Release Retractor Handle	2
SB-2010	Clear-Line Mini Sterilization Container	1
SB-2330B	30mm Blunt Blade	2
SB-2340B	40mm Blunt Blade	2
SB-2350B	50mm Blunt Blade	2
SB-2360B	60mm Blunt Blade	2
SB-2370B	70mm Blunt Blade	2
SB-2330N	30mm Notched Blade	2
SB-2340N	40mm Notched Blade	2
SB-2350N	50mm Notched Blade	2
SB-2360N	60mm Notched Blade	2
SB-2370N	70mm Notched Blade	2

 V797D-70161



Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support

Main +1 (800) 303-3050 Fax +1 (949) 600-8874 9620 Research Drive, Irvine - CA 92618 USA info@sbhsurgical.com www.sbhsurgical.com

Copyright © SBH ®. 2021 - All rights reserved.