



Strategic Business Holdings, Inc.


...Quality Products, Superior Customer Service, & Support




SBH 3380 Fechtner Conjunctiva Forceps
very delicate ring shaped tips 1.2mm
for optimal: grasping without Cutting or Tearing



SBH 6093 Mermound Forceps
very delicate tips
for optimal: micro manipulation



SBH 6094 Inamura Capsulorrhexis Forceps
cross-action
for optimal: vision and controlled action



SBH 6082 Crozafon Capsulorrhexis Forceps
cystotome shaped tip, maximum jaw opening 2 mm
optimal for: for 1.5 mm incisions



SBH 6083 Gianetti Capsulorrhexis Forceps
vaulted thin shanks
for optimal: access without stretching the incision



SBH 6082 Ernest Nucleus Cracker
cross action, delicate, serrated paddle tips
for optimal: dividing of nucleus into smaller segments

100% made in Germany

Capsulorrhexis Forceps