

The GRAY Laparoscopic Upper GI Assistant ensures quick, stable, unobtrusive scope and instrument fixation for operations on the upper stomach and is specifically designed for:

- Gastric Banding and
- Nissen Fundoplication



## Laparoscopic Upper GI Assistant

1. MH01T Multiflex Holder Teflon
2. LL 30 / LL50 / LL70 Liver Lift Set
3. BC01 Coupling Clamp
4. LHA01 Extension Arm
5. EC05 Grip Clamp
6. GT01 Table Clamp



GRAY SURGICAL

4



The **GRAY Laparoscopic Upper GI Assistant** is extremely user friendly because it:

- (a) elevates the left lobe of the liver out of the surgical field and
- (b) holds the endoscope/light source in a stable position while allowing it to be instantly moved to any desired position and
- (c) removes the need for an assistant during these operations



Quick and easy to set up



Teflon jaws allow ease of scope manipulation

The **GRAY Laparoscopic Upper GI Assistant** is designed to ensure flexible, multidirectional scope adjustment and effective liver retraction.

- \* Holds the endoscope/light source securely in position
- \* Eliminates hand held retraction
- \* Minimises hand held scope manipulation
- \* Teflon jaws allow ease of scope manipulation and eliminates risk of damage to scope shaft
- \* A set of 3 uniquely designed Liver Lift retractors to retract the liver and provide even distribution of the weight load
- \* Quick and easy to set up
- \* GT01 Table clamp fits all operating table rails
- \* Warranted for one year against manufacturing defects



LL30/ LL50/ LL70 Liver Lift Set

Patents Apply