

The GRAY Endoscopic Assistant ensures quick, stable, unobtrusive scope and instrument fixation for all endoscopic surgical procedures.

The Gray Endoscopic Assistant is provided with either **Delrin Jaws** for holding instruments or **Teflon Jaws** for holding endoscopes/light sources.



2

1



**Single Arm System**

- 1. MH01T Multiflex Holder Teflon \*
- 2. GT01 Table Clamp

For more information you may contact us at:



**Strategic Business Holdings, Inc.**

9620 Research Dr. Irvine, CA 92618 USA

Tel. 800-303-3050 Fax 949-600-8847

[www.sbhsurgical.com](http://www.sbhsurgical.com) [info@sbhsurgical.com](mailto:info@sbhsurgical.com)



Multiple 'Assistants' can be used during the one surgical procedure and can be mounted either on the Table Post or Extension Arm.



Delrin jaws - recommended for use with instruments



Teflon jaws - recommended for use with endoscopes/light sources

### The GRAY Endoscopic Assistant is designed to ensure flexible, multidirectional adjustment and fixation.

- \* Holds any instrument, retractor or endoscope securely in position
- \* Eliminates hand held retraction
- \* Minimises hand held scope manipulation
- \* Provides a secure grip on 3mm-12mm instrument shafts
- \* Teflon jaws allow ease of scope manipulation and eliminates risk of damage to scope shaft
- \* Delrin jaws provide secure grip of laparoscopic instruments
- \* Quick and easy to set up
- \* GT01 Table clamp fits all operating table rails
- \* Warranted for one year against manufacturing defects

### **\*Important:**

When ordering the Multiflex Holder please specify either:  
 a) MH01T Teflon - recommended for use with endoscopes/light sources  
 b) MH01D Delrin - recommended for use with instruments

Patents Apply