

The GRAY Laparoscopic Cholecystectomy Assistant ensures quick, stable, unobtrusive scope and instrument fixation during laparoscopic cholecystectomy.



The **GRAY Laparoscopic Cholecystectomy Assistant** is extremely user friendly because it:

- (a) elevates the gallbladder to expose the operative field and
- (b) holds the endoscope/light source in a stable position while allowing it to be instantly moved to any desired position and
- (c) removes the need for an assistant during these operations



**Laparoscopic Cholecystectomy Assistant**

- 1. MH01D Multiflex Holder Delrin
- 2. MH01T Multiflex Holder Teflon
- 3. CC01 Coupling Clamp
- 4. GE01 Extension Arm
- 5. GT01 Table Clamp



GRAY SURGICAL



**The Gray Laparoscopic Cholecystectomy Assistant is designed to ensure flexible, multidirectional adjustment and fixation.**

- \* Holds any instrument, retractor or endoscope securely in position
- \* Eliminates hand held retraction
- \* Minimises hand held scope manipulation
- \* Provides a secure grip on 3mm-12mm instrument shafts
- \* Teflon jaws allow ease of scope manipulation and eliminates risk of damage to scope shaft
- \* Delrin jaws provide secure grip of laparoscopic instruments
- \* Quick and easy to set up
- \* GT01 Table clamp fits all operating table rails
- \* Warranted for one year against manufacturing defects

**\*Important:**

When ordering the Multiflex Holder please specify either:  
a) MH01T Teflon - recommended for use with endoscopes/light sources  
b) MH01D Delrin - recommended for use with instruments



Quick and easy to set up



Delrin jaws - recommended for use with instruments



Teflon jaws - recommended for use with endoscopes / light sources

Patents Apply