




Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support





100% made in Germany

MF Cranial Stabilization Systems

Copyright © SBH ®. 2026 - All rights reserved.



MF Cranial Stabilization Systems





SBH: Your address for surgical quality instruments "Made in Germany", manufactured at our site in Germany. SBH stands for highest quality, speed, and efficiency. We will find the technically and economically optimal solution for you.

Highest quality at fair prices. Contact us to get to know the complete SBH product range!

For further information please visit: www.sbhsurgical.com for by e-mail at info@sbhsurgical.com

1 MF- Functional overview	4-5
2 MF-Headrest System.....	6
3 MF- Skull Pins.....	7
4 MF-Horseshoe Headrest.....	8-9
5 MF-Retractor System	10-11
6 Table Adapters.....	12-13
7 MF-Retractor System (Child Size)	14
Contact	16

MF-Cranial Stabilization System

Material

- MF-Cranial fixation system made of aluminium I stainless steel
- Coated MF-Cranial fixation system
- MF-Cranial fixation system made of high-tech plastic for complete X-ray transparency and MRI-compatibility

Product Design

- The gold standard in head fixation
- 2 sets available with horseshoe head-rest (for adults) or with skull clamp (3-pin-jixation)
- MF-Pins in stainless steel, titanium or PEEK

Description

- As a complete system with all components, Made in Germany
- Outstanding stability with all necessary flexibility
- Components made of aluminium for weight reduction or stainless steel for stability
- Milled components for improved stability
- Basic components compatible with existing systems
- Can be combined with Viridis and snake-arm-system




MF-Cranial Stabilization System

Functional Overview



N = Navigation



Availability

[AL] Gold Standard

[AL] + [SST] Premium coated

[CT] CT compatible

[MRI] MRI compatible

Unique solution

The newly developed MF head stabilization and retraction system has been specially designed to meet the needs of the user.

The basic set is available with a skull clamp or neck support and can be connected to a complete MF attachment holding system for the adaptation of accessories (Viridis, Leyla ..).

This makes the system movable in all directions and height-adjustable. Four adaption interfaces for the attachment of target devices as well as further connection possibilities for stereotactic procedures are available.

Everything from a single source. Made in Germany with the highest quality.

Can be adapted with two bars per side for more flexibility. Can also be used as a hand rest.

Variants: straight or curved rods in various lengths.

Round, narrow contour of the rods: can be adapted with standards such as Leyla/Viridis.

Navigation adapter connection.

Protection against unintentional release of the locking device.

Four adaptation points available for target devices / navigation adapter connection (N). Further connection options for stereotactic procedures.

Individual labelling is possible here.

Clamping adapter 360° rotatable and tiltable.

Production from full aluminium block for higher stability and smoother pore-free surface.

Operating elements color marked (black/gold).



MF-Headrest System



Aluminium

Table adapters
on pages 16-17



SBH 117-000-10

MF-Headrest System with skull clamp

consisting of:

SBH 117-000-15

SBH 117-000-25

SBH 117-000-40



SBH 117-000-40

MF-Skull Clamp



SBH 117-000-15

MF-Adjustable Base Unit with connecting arm



SBH 117-000-25

MF-Swivel Adaptor

MF-Skull Pins



Stainless Steel **SST**

Titanium **TI**

Peek **PE**



Article Number	Description	Material		Application	Adult/Pediatric
SBH 117-000-50	MF-Skull Pins	Stainless Steel	SST	Reusable	Adults
SBH 117-000-52	MF-Skull Pins	Stainless Steel	SST	Reusable	Pediatric
SBH 117-000-54	MF-Skull Pins	Titanium	TI	Reusable	Adults
SBH 117-000-56	MF-Skull Pins	Titanium	TI	Reusable	Pediatric
SBH 117-000-62	MF-Skull Pins	PEEK + Stainless Steel	PE SST	Disposable	Adults
SBH 117-000-64	MF-Skull Pins	PEEK + Stainless Steel	PE SST	Disposable	Pediatric
SBH 117-000-66	MF-Skull Pins	PEEK + Titanium	PE TI	Disposable	Adults
SBH 117-000-68	MF-Skull Pins	PEEK + Titanium	PE TI	Disposable	Pediatric

Packaging unit: 1 pc



SBH 117-000-50



SBH 117-000-52



SBH 117-000-54



SBH 117-000-56



SBH 117-000-62




SBH 117-000-64



SBH 117-000-66



SBH 117-000-68





MF-Horseshoe Headrest



Aluminium



SBH 117-000-05

MF-Headrest System with Horseshoe Headrest

consisting of:

SBH 117-000-15

SBH 117-000-25

SBH 117-000-35

Aluminium



Table adapters
on page 16-17



SBH 117-000-35

MF-Horseshoe Headrest (incl. gel pad)



SBH 117-000-37

MF-Gel Pad Adults (set)



SBH 117-000-38

MF-Horseshoe Headrest Basis (without gel pad)




SBH 117-000-15

MF-Adjustable Base Unit with connecting arm



SBH 117-000-25

MF-Swivel Adaptor






MF-Retractor System

F1

Stainless steel



SBH 117-000-70	MF-Retractor System for adaption of accessories (Viridis, Leyla...) consisting of:	
Set content	Description	Quantity
SBH 117-000-71	MF-Connection arm to basic headrest system (Set, right-left)	1
SBH 116-015-00	MF-locking connector for clamp and rod	2
SBH 117-000-73	MF-support rod with gear	2
SBH 117-000-74	MF-arm Ø12,7 slightly curved (Set, right-left)	2
SBH 117-000-72	MF-Retractor Set-up (1 side), right-left adaptable consisting of:	
Set content	Description	Quantity
SBH 117-000-71	MF-Connection arm to basic headrest system (Set, right-left)	1
SBH 116-015-00	MF-locking connector for clamp and rod	1
SBH 117-000-73	MF-support rod with gear	1
SBH 117-000-74	MF-arm Ø 12, 7 slightly curved (Set, right-left)	1

MF-Retractor System



Stainless steel

SST

Aluminium

AL



SBH 117-000-71

MF-Connection arm to basic headrest system (Set, right-left)




SBH 116-015-00

MF-locking connector for clamp and rod



SBH 117-000-74

MF-arm Ø 12,7 slightly curved (Set, right and left)



SBH 117-000-76

MF-arm Ø 12,7, short, straight (Set, right and left)

SBH 117-000-75

MF-arm Ø 12,7 curved (Set, right and left)

SBH 117-000-77

MF-arm Ø 12,7, long, straight (Set, right and left)



SBH 117-000-73

MF-support rod with gear



SBH 117-001-35

MF-Universal Slide Rail Fitting
(for MB1-17-001-00, please order 2 psc.)



SBH 117-001-00

MF-Crossbar, width adaptable (64 - 70 cm)






Table Adapters

Combination options and selection of table adapters



Table adapter with round tube and short side arms with rails

SBH 117-000-32 Model S

SBH 117-000-42 Model M

SBH 117-000-44 Model L



Table adapter with square rod and short side arms with rails

SBH 117-000-31 Model S

SBH 117-000-33 Model M

SBH 117-000-43 Model L




Table adapter Porkbench with long short side arms and round tube with rails

SBH 117-001-20 Model S


SBH 117-001-21 Model M

SBH 117-001-22 Model L




SBH 117-002-90

Connection spigot 025.4mm (1")
(pair) + 2 x screw for e.g. Trump card



SBH 117-002-92

Connection spigot 019 (pair) + 2 x screw



SBH 117-002-94

Connection spigot 0.18 (pair) + 2 x screw
for e.g. Moquet Yuno II



SBH 117-001-40

MF-Armrest Porkbench



Table adapters for MF-headrest system and premium MF-headrest system



Article number	Description	Rail distance in mm	Bore spacing in mm
SBH 117-000-30	MF-Siderail OR table adapter	550-610	



MF-Table adapter universal with round tube and short side arms with rails

Article Number	Description	Rail Distance in mm (R)	Bore Spacing in mm (B)
SBH 117-000-32	Model S (small size)	366-531	386-551
SBH 117-000-42	Model M (middle size)	526-690	546-710
SBH 117-000-44	Model L (large size)	685-850	705-870

To be used with **SBH 117-100-20/25/40/35**


MF- Side Rail OR Table adapter universal with square rod with rails.

Article Number	Description	Rail distance in mm (R)	
SBH 117-000-31	Model S (small size)	400-520	420-540
SBH 117-000-33	Model M (middle size)	530-650	550-670
SBH 117-000-43	Model L (large size)	730-850	750-870

To be used with **SBH 117-000-15/25/40/35** or **SBH 117-100-15/25/40/35**

Table adapter Parkbench with long short side arms and around tube with rails.


Article number	Description	Rail distance in mm (R)	Bore spacing in mm (B)
SBH 117-001-20	MF-Adjustable Base Unit Parkbench, small S	366-531	386-551
SBH 117-001-21	MF-Adjustable Base Unit Parkbench, medium M	526-690	546-710
SBH 117-001-22	MF-Adjustable Base Unit Parkbench, long L	685-850	705-870

To be used with **SBH 117-000-15/25/40/35** or **SBH 117-100-15/25/40/35**





MF Cranial Stabilization System

MF Cranial Stabilization Systems Child Size



SBH 117-200-45



[illegible]



V797D-70161



Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support

Main +1 (800) 303-3050 Fax +1 (949) 600-8874 9620 Research Drive, Irvine - CA 92618 USA info@sbhsurgical.com www.sbhsurgical.com

Copyright © SBH ®. 2026 - All rights reserved.