




Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support



Sellors Rib Spreader

- Curved arms
- Keyless adaption of blades to frame
- Stainless steel
- 100% made in Germany

SELLORS RIB SPREADER COMPLETE

SB 110-739-05	consisting of:			
---------------	----------------	--	--	--

	Type	spreading width	arm length	Qty
SB 110-739-10	Sellors Frame	min. 197mm max. 226mm	140mm	1

	Type	width	depth	Qty
SB 110-739-25	1 pair of blades	50 mm	40 mm	1
SB 110-739-30	1 pair of blades	50 mm	50 mm	1
SB 110-739-35	1 pair of blades	50 mm	70 mm	1

OPTIONS

	Type	width	depth	Qty
SB 110-739-50	1 pair of blades	75 mm	30 mm	1
SB 110-739-55	1 pair of blades	75 mm	45 mm	1
SB 110-739-60	1 pair of blades	120 mm	30 mm	1
SB 110-739-65	1 pair of blades	75 mm	78 mm	1

100% made in Germany

Sellors Rib Spreader

Main +1 (800) 303-3050 Fax +1 (949) 600-8874 9620 Research Drive, Irvine - CA 92618 USA info@sbhsurgical.com www.sbhsurgical.com

Copyright © SBH ©. 2020 - All rights reserved.