




Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support

- Lumbar retractor system
- Side load for optimal exposure
- Titanium blades, color-coded, radiolucent
- Large variety of valves and hooks




100% made in Germany

CLR Retractor Set



CLR Retractor Set



SBH 110-751-05*

CLR retractor side-load



SBH 110-804-00*

CCR/ CLR handle for blades with black PPSU-handle



SBH 110-804-05*

CCR/ CLR handle for blades, cannulated with Luer-Lock, with black PPSU-handle




SBH 110-794-10*

CCR / CLR depth gauge instrument, natural

CLR Retractor Set



All blades out of Titanium, color: blue, radiolucent, with side-load mechanism



CLR blade, side load, width: 20mm, blue

	B - width	T - depth
SBH 110-749-21*	20 mm	35 mm
SBH 110-749-23*	20 mm	40 mm
SBH 110-749-25*	20 mm	45 mm
SBH 110-749-27*	20 mm	55 mm
SBH 110-749-29*	20mm	65 mm




CLR blade, side load, width: 40mm, blue

	B - width	T - depth
SBH 110-749-31*	40 mm	35 mm
SBH 110-749-33*	40 mm	40 mm
SBH 110-749-35*	40 mm	45 mm
SBH 110-749-37*	40 mm	55 mm
SBH 110-749-39*	40 mm	65 mm
SBH 110-749-41*	40 mm	75 mm




CLR blade, side load, width: 60mm, blue


	B - width	T - depth
SBH 110-749-43*	60 mm	35 mm
SBH 110-749-45*	60 mm	40 mm
SBH 110-749-47*	60 mm	45 mm
SBH 110-749-49*	60 mm	55 mm
SBH 110-749-51*	60 mm	65 mm
SBH 110-749-53*	60 mm	75 mm
SBH 110-749-55*	60 mm	85 mm



Hemi-laminectomy hook, blue


SBH 110-749-61*	small
SBH 110-749-66*	large





SBH 120-968-73*

CLR-Tray, stainless steel with cover, dimensions: 485 x 253 x 70 mm



SBH 110-751-10

CLR Set side load consisting of products marked (*)
Blades, 2 pcs each

Changes in technical specifications and dimensions are at our discretion. Subject to modification.
Without our explicit written approval this print media may not be copied or reprinted as a whole or in extracts.



Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support

Main +1 (800) 303-3050 Fax +1 (949) 600-8874 9620 Research Drive, Irvine, CA 92618 USA info@sbhsurgical.com www.sbhsurgical.com

Copyright © SBH ©. 2020 - All rights reserved.