




# Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support



## KC Frame Retractor Set

- For posterior lumbar and thoracic spine surgery.
- Standard Set with KC Frame Retractor and standard blades in different designs.
- 1 to 5 teeth design, various widths and depths of blades available.
- Optimum assembly, repositioning and stability.
- All instruments available in one tray, securely stored (dimensions: 540 x 253 x 50mm).

## STANDARD KC FRAME RETRACTOR SET

<b>SBH 110-762-05</b>	without tray, consisting of:			
<b>SBH 110-762-06</b>	with tray and lid <b>SBH 120-968-78</b> , consisting of:			
	<i>type</i>	<i>width</i>	<i>length</i>	<i>quantity</i>
<b>SBH 110-762-10</b>	<b>Standard KC Frame Retractor Set</b>			1
<b>SBH 110-764-15</b>	5 teeth	76 mm	64 mm	2
<b>SBH 110-764-03</b>	left	31 mm	38 mm	1
<b>SBH 110-764-06</b>	right	31 mm	38 mm	1
<b>SBH 110-764-09</b>	left	57 mm	44 mm	1
<b>SBH 110-764-12</b>	right	57 mm	44 mm	1

*KC-Blades standard*



100% made in Germany

## KC Frame Retractor Set With Tray



# Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support



New:

- **Lightweight K Frame Retractor Set with newly designed retractor** and connection of blades with additional safeguard.
- Retractor is 25% lighter compared to standard retractor due to new design.
- Blades feature an additional safeguard to allow fixation of blades.

## LIGHTWEIGHT KC FRAME RETRACTOR SET

<b>SBH 110-762-13</b>	without tray, consisting of:
<b>SBH 110-762-14</b>	with tray and lid <b>SBH 120-968-78</b> , consisting of:

	<i>type</i>	<i>width</i>	<i>length</i>	<i>quantity</i>
<b>SBH 110-762-11</b>	<b>Lightweight KC Frame Retractor Set</b>			1
<b>SBH 110-764-65</b>	5 teeth	76 mm	64 mm	2
<b>SBH 110-764-54</b>	left	31 mm	38 mm	1
<b>SBH 110-764-56</b>	right	31 mm	38 mm	1
<b>SBH 110-764-59</b>	left	57 mm	44 mm	1
<b>SBH 110-764-62</b>	right	57 mm	44 mm	1

*KC-Blades with special safeguard*

<b>SBH 120-968-78</b>	tray with lid, 540 x 253 x 50 mm
-----------------------	----------------------------------

100% made in Germany

## KC Frame Retractor Set With Tray