

# LIVER RETRACTOR

The exposure presented by this unique design is unrivalled in Liver surgery.



## GRAY SURGICAL

Precision Is Our Business

Manufacturing Quality Retraction  
Systems For The Specialist Surgeon

[www.sbhsurgical.com](http://www.sbhsurgical.com)

JAS-ANZ



CE

TGA

Patents Apply

Issue - 0102

Code - LRB1

## Quality

- Precision made
- Highest quality
- Superior warranty

## Versatility

- One design fits all operating tables.
- Superior stability.
- Ease of assembly
- Fully adjustable Upper Hoop.
- Unique self-aligning blades allows for maximum rib retraction presenting the liver centrally within the surgical field.
- Excellent surgical exposure
- Quick, secure, blade attachment.
- Numerous blade designs offer versatility of retraction.
- Customised blade designs available

## Convenience

- Eliminates the need for manual retraction
- Multiple Blade configurations


### Components




Liver Table Clamp  
-code LT01 (2 req)  
inc 10 Teflon Insulators  
-code TE01



Upper Hoop  
-code UH01



Risers  
-code RI01 (2 req)



Lower Hoop  
-code LH01



Liver Wedges  
-code WD01 (6 req)

Blades - complete range available - refer to Blades Brochure

For complete range of Gray Surgical Products, please see our website at [www.graysurgical.com](http://www.graysurgical.com)

Patents

Apply



Your Local Distributor:

**Strategic Business Holdings, Inc.**  
...quality products, superior customer service & support

Tel: 800.303.3050 Fax: 949.600.8874



GRAY SURGICAL

Please refer to our website for Warranty and Care and Maintenance guidelines of Gray Surgical Equipment.  
For more information on Gray Surgical Retractor Systems please contact your local distributor or visit our website

9620 Research Dr. Irvine, Ca 92618 USA

[www.sbhsurgical.com](http://www.sbhsurgical.com)

[info@sbhsurgical.com](mailto:info@sbhsurgical.com)