




Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support

Micro Discectomy Retractor




- Optimal Exposure
- Retractors and blades with TiN coated black finish

100% made in Germany

MCC Retractor System



MCC Retractor System



Retractors: TiN coated

SB 110-777-05	Set SB 110-777-20	Set SB 110-777-25
MCC Hinged Retractor max. opening: 85 mm/ 3.3"	*	**

SB 110-777-10

**MCC Rigid Retractor
max. opening: 85 mm/ 3.3"**

All blades and hooks: TiN coated

Blade 2 teeth, 17 mm

	B - width	T - depth
SB 110-778-82	17 mm	30 mm
SB 110-778-84	17 mm	40 mm
SB 110-778-86	17 mm	50 mm
SB 110-778-88	17mm	60 mm
SB 110-778-90	17 mm	70 mm
SB 110-778-92	17 mm	80 mm
SB 110-778-93	17 mm	90 mm

Blade 2 teeth, 20 mm

	B - width	T - depth
SB 110-778-42	20 mm	50 mm
SB 110-778-43	20 mm	60 mm
SB 110-778-44	20 mm	70 mm

Tapered blade, 6 teeth, 27 mm

	B - width	T - depth
SB 110-778-95	27 mm	40 mm

MCC Retractor System



Blade, 4 teeth, 20 mm

	B - width	T - depth	Set SB 110-777-20	Set SB 110-777-25
SB 110-778-05	20 mm	30 mm		**
SB 110-778-10	20 mm	40 mm	*	**
SB 110-778-15	20 mm	50 mm	*	**
SB 110-778-20	20 mm	60 mm	*	**
SB 110-778-25	20 mm	70 mm	*	**
SB 110-778-30	20 mm	80 mm		**
SB 110-778-35	20 mm	90 mm		



Blade, 6 teeth, 27 mm

	B - width	T - depth	Set SB 110-777-20	Set SB 110-777-25
SB 110-778-55	27 mm	30 mm		**
SB 110-778-60	27 mm	40 mm	*	**
SB 110-778-65	27 mm	50 mm	*	**
SB 110-778-70	27 mm	60 mm	*	**
SB 110-778-75	27 mm	70 mm	*	**
SB 110-778-80	27 mm	80 mm		**
SB 110-778-81	27 mm	90 mm		




Multiple toothed hook , 2.0 mm

	Length	Tooth	Set SB 110-777-20	Set SB 110-777-25
SB 110-779-55	20 mm	2		**
SB 110-779-60	30 mm	3	*	**
SB 110-779-65	40 mm	4	*	**
SB 110-779-70	50 mm	5	*	**
SB 110-779-75	60 mm	6	*	**
SB 110-779-80	70 mm	7		**
SB 110-779-81	80 mm	8		



Hook blade, 1 teeth, 1.8 mm

	Length	Set SB 110-777-20	Set SB 110-777-25
SB 110-779-05	20 mm		**
SB 110-779-10	30 mm	*	**
SB 110-779-15	40 mm	*	**
SB 110-779-20	50 mm	*	**
SB 110-779-25	60 mm	*	**
SB 110-779-30	70 mm		**
SB 110-779-35	80 mm		





MCC Retractor System

SB 120-968-74 *


MCC Tray small with lid 243 mm x 253 mm x 50 mm

SB 120-968-75 **

MCC Tray large with lid 405 mm x 253 mm x 50 mm


SB 110-777-20

Basic MCC Set with tray and hinged retractor (*)
Blades 2 each, hook 1 each



SB 110-777-25

MCC Set large with tray
and 2 hinged retractors (**)
Blades 2 each, hook 1 each



Changes in technical specifications and dimensions are at our discretion. Subject to modification.
Without our explicit written approval this print media may not be copied or reprinted as a whole or in extracts.



Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support

Main +1 (800) 303-3050 Fax +1 (949) 600-8874 9620 Research Drive, Irvine - CA 92618 USA info@sbhsurgical.com www.sbhsurgical.com

Copyright © SBH ©. 2020 - All rights reserved.