



Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support


Sterilization Container Systems



- Safe, reliable, easy to handle
- Lightweight
- Energy saving



Environmental
Friendly





Sterilization Container Systems

Page

<u>02</u>	<u>Index</u>
<u>03-04</u>	<u>Information</u>
<u>05-11</u>	<u>1/1 Container and Wire Baskets</u>
<u>12-17</u>	<u>3/4 Container and Wire Baskets</u>
<u>18-25</u>	<u>1/2 Container and Wire Baskets</u>
<u>26-27</u>	<u>Mini Container</u>
<u>28-29</u>	<u>Dental Container</u>
<u>31-32</u>	<u>Accessories</u>

NOTES:

All our sterilization containers have the sterilization validation documentation that are approved by the FDA and have a FDA Class II 510 (k).

Standard Model

Choose between our Standard and our Safety Model.

The safety lid protects against contamination during transportation and storage.

The filter and standard safety lid can only be separated by complete opening of the container.

For a simple detachment of the safety lid from the perforated lid we offer also a new model of safety lid with exterior press-button.



New Safety Model



Sterilization Container Systems

Information

Benefits

The new container models with filter are provided with a longer handle, which guarantees an easy and safe manipulation. The silicone tube around the handle protects against burning and sliding.



The closures of these containers are furnished with a feather system that insures a stable closing and cancels the overpressure that is formed within the container.



The filter retainer is locked by pressing down in center of the filter. The opening of the filter retainer succeeds by pressing the plate.



Container Size

111



Standard Model

New Safety Model



Sterilization Container Systems

Standard Model

Please note that only the lid of the container is coloured.

1/1



100

outside
580 x 280 x 100 mm
inside
550 x 265 x 84 mm

Bottom non-perforated



Bottom perforated



SB 110.100E
SB 120.100E
SB 130.100E
SB 140.100E
SB 150.100E
SB 160.100E

SB 111.100E
SB 121.100E
SB 131.100E
SB 141.100E
SB 151.100E
SB 161.100E



135

outside
580 x 280 x 135 mm
inside
550 x 265 x 115 mm

SB 110.135E
SB 120.135E
SB 130.135E
SB 140.135E
SB 150.135E
SB 160.135E

SB 111.135E
SB 121.135E
SB 131.135E
SB 141.135E
SB 151.135E
SB 161.135E



150

outside
580 x 280 x 150 mm
inside
550 x 265 x 130 mm

SB 110.150E
SB 120.150E
SB 130.150E
SB 140.150E
SB 150.150E
SB 160.150E

SB 111.150E
SB 121.150E
SB 131.150E
SB 141.150E
SB 151.150E
SB 161.150E



200

outside
580 x 280 x 200 mm
inside
550 x 265 x 180 mm

SB 110.200E
SB 120.200E
SB 130.200E
SB 140.200E
SB 150.200E
SB 160.200E

SB 111.200E
SB 121.200E
SB 131.200E
SB 141.200E
SB 151.200E
SB 161.200E




260

outside
580 x 280 x 260 mm
inside
550 x 265 x 240 mm

SB 110.260E
SB 120.260E
SB 130.260E
SB 140.260E
SB 150.260E
SB 160.260E

SB 111.260E
SB 121.260E
SB 131.260E
SB 141.260E
SB 151.260E
SB 161.200E

Bottom non-perforated



Bottom perforated



1/1

Standard Safety Model

Please note that only the lid of the container is coloured.

SB 112.100E
SB 122.100E
SB 132.100E
SB 142.100E
SB 152.100E
SB 162.100E

SB 113.100E
SB 123.100E
SB 133.100E
SB 143.100E
SB 153.100E
SB 163.100E

100

outside
580 x 280 x 100 mm
inside
550 x 265 x 84 mm



SB 112.135E
SB 122.135E
SB 132.135E
SB 142.135E
SB 152.135E
SB 162.135E

SB 113.135E
SB 123.135E
SB 133.135E
SB 143.135E
SB 153.135E
SB 163.135E

135

outside
580 x 280 x 135 mm
inside
550 x 265 x 115 mm



SB 112.150E
SB 122.150E
SB 132.150E
SB 142.150E
SB 152.150E
SB 162.150E

SB 113.150E
SB 123.150E
SB 133.150E
SB 143.150E
SB 153.150E
SB 163.150E

150

outside
580 x 280 x 150 mm
inside
550 x 265 x 130 mm



SB 112.200E
SB 122.200E
SB 132.200E
SB 142.200E
SB 152.200E
SB 162.200E

SB 113.200E
SB 123.200E
SB 133.200E
SB 143.200E
SB 153.200E
SB 163.200E

200

outside
580 x 280 x 200 mm
inside
550 x 265 x 180 mm



SB 112.260E
SB 122.260E
SB 132.260E
SB 142.260E
SB 152.260E
SB 162.260E

SB 113.260E
SB 123.260E
SB 133.260E
SB 143.260E
SB 153.260E
SB 163.260E

260

outside
580 x 280 x 260 mm
inside
550 x 265 x 240 mm





Sterilization Container Systems

New Safety Model

Variant with Quick Opening Mechanism

Please note that only the lid of the container is coloured.

1/1

Bottom non-perforated



Bottom perforated



100

outside
580 x 280 x 100 mm
inside
550 x 265 x 84 mm

SB 116.100E
SB 126.100E
SB 136.100E
SB 146.100E
SB 156.100E
SB 166.100E

SB 117.100E
SB 127.100E
SB 137.100E
SB 147.100E
SB 157.100E
SB 167.100E



135

outside
580 x 280 x 135 mm
inside
550 x 265 x 115 mm

SB 116.135E
SB 126.135E
SB 136.135E
SB 146.135E
SB 156.135E
SB 166.135E

SB 117.135E
SB 127.135E
SB 137.135E
SB 147.135E
SB 157.135E
SB 167.135E



150

outside
580 x 280 x 150 mm
inside
550 x 265 x 130 mm

SB 116.150E
SB 126.150E
SB 136.150E
SB 146.150E
SB 156.150E
SB 166.150E

SB 117.150E
SB 127.150E
SB 137.150E
SB 147.150E
SB 157.150E
SB 167.150E



200

outside
580 x 280 x 200 mm
inside
550 x 265 x 180 mm

SB 116.200E
SB 126.200E
SB 136.200E
SB 146.200E
SB 156.200E
SB 166.200E

SB 117.200E
SB 127.200E
SB 137.200E
SB 147.200E
SB 157.200E
SB 167.200E



260

outside
580 x 280 x 260 mm
inside
550 x 265 x 240 mm

SB 116.260E
SB 126.260E
SB 136.260E
SB 146.260E
SB 156.260E
SB 166.260E

SB 117.260E
SB 127.260E
SB 137.260E
SB 147.260E
SB 157.260E
SB 167.200E

1/1

Bottoms and Lids

Bottom non perforated




Article-No.	outside	inside
SB 105.100E	580 x 280 x 100 mm	550 x 262 x 84 mm
SB 105.135E	580 x 280 x 135 mm	550 x 262 x 117 mm
SB 105.150E	580 x 280 x 150 mm	550 x 262 x 132 mm
SB 105.200E	580 x 280 x 200 mm	550 x 262 x 182 mm
SB 105.260E	580 x 280 x 260 mm	550 x 262 x 244 mm

Bottom perforated



Article-No.	outside	inside
SB 106.100E	580 x 280 x 100 mm	550 x 262 x 84 mm
SB 106.135E	580 x 280 x 135 mm	550 x 262 x 117 mm
SB 106.150E	580 x 280 x 150 mm	550 x 262 x 132 mm
SB 106.200E	580 x 280 x 200 mm	550 x 262 x 182 mm
SB 106.260E	580 x 280 x 260 mm	550 x 262 x 244 mm

silver
yellow
red
blue
green
black



580 x 280 mm

Lid perforated




SB 100.010E
SB 100.020E
SB 100.030E
SB 100.040E
SB 100.050E
SB 100.060E

Safety Lid- New Model- with exterior press- Button



580 x 280 mm

Safety Lid



Standard	Standard
SB 103.010	SB 103.010.E
SB 103.020	SB 103.020.E
SB 103.030	SB 103.030.E
SB 103.040	SB 103.040.E
SB 103.050	SB 103.050.E
SB 103.060	SB 103.060.E



Accessories 1/1




Stainless steel wire basket
intermediate crimp square hole wire mesh

Article-No.	Dimensions
SB 191.030	540 x 255 x 35 mm
SB 191.050	540 x 255x 50 mm
SB 191.070	540 x 255 x 70 mm
SB 191.100	540 x 255 x 100 mm




Stainless steel wire basket
intermediate crimp square hole wire mesh

Article-No.	Dimensions
SB 191.030	485 x 255 x 30 mm
SB 191.050	485 x 255 x 50 mm
SB 191.070	485 x 255 x 70 mm
SB 191.100	485 x 255 x 100 mm




Stainless steel screen basket
with lid

Article-No.	Dimensions
SBS 191.030.02	545 x 255 x 35 mm
SBS 191.050.02	545 x 255 x 55 mm
SBS 191.070.02	545 x 255 x 75 mm
SBS 191.100.02	545 x 255 x 105 mm



Lid for stainless steel screen basket
with closure

Article-No.	Dimensions
SBS 191.000	545 x 255 mm



Stainless steel screen basket

Article-No.	Dimensions
SBS 191.030	545 x 255 x 30 mm
SBS 191.050	545 x 255 x 50 mm
SBS 191.070	545 x 255 x 70 mm
SBS 191.100	545 x 255 x 100 mm

1/1

Accessories

Article-No.	Dimensions
SBS 190.030.02	485 x 255 x 35 mm
SBS 190.050.02	485 x 255 x 55 mm
SBS 190.070.02	485 x 255 x 75 mm
SBS 190.100.02	485 x 255 x 105 mm

Stainless steel screen basket with lid



Article-No.	Dimensions
SBS 190.000	485 x 255 mm

Lid for stainless steel screen basket with closure



Article-No.	Dimensions
SBS 190.030	485 x 255 x 30 mm
SBS 190.050	485 x 255 x 50 mm
SBS 190.070	485 x 255 x 70 mm
SBS 190.100	485 x 255 x 100 mm

Stainless steel screen basket





Container Size


314



Standard Model

New Safety Model

Bottom non-perforated



Bottom perforated



SB 210.100E
SB 220.100E
SB 230.100E
SB 240.100E
SB 250.100E
SB 260.100E

SB 211.100E
SB 221.100E
SB 231.100E
SB 241.100E
SB 251.100E
SB 261.100E

100

outside
465 x 280 x 100 mm
inside
425 x 260 x 84 mm



SB 210.135E
SB 220.135E
SB 230.135E
SB 240.135E
SB 250.135E
SB 260.135E

SB 211.135E
SB 221.135E
SB 231.135E
SB 241.135E
SB 251.135E
SB 261.135E

135

outside
465 x 280 x 135 mm
inside
425 x 260 x 115 mm



SB 210.150E
SB 220.150E
SB 230.150E
SB 240.150E
SB 250.150E
SB 260.150E

SB 211.150E
SB 221.150E
SB 231.150E
SB 241.150E
SB 251.150E
SB 261.150E

150

outside
465 x 280 x 150 mm
inside
425 x 260 x 130 mm



3/4

Standard Model

Please note that only the lid of the container is coloured.



Sterilization Container Systems

Standard Safety Model

Please note that only the lid of the container is coloured.


3/4



100

outside
465 x 280 x 100 mm
inside
425 x 260 x 84 mm

Bottom non-perforated



Bottom perforated



SB 212.100E
SB 222.100E
SB 232.100E
SB 242.100E
SB 252.100E
SB 262.100E

SB 213.100E
SB 223.100E
SB 233.100E
SB 243.100E
SB 253.100E
SB 263.100E



135

outside
465 x 280 x 135 mm
inside
425 x 260 x 115 mm

SB 212.135E
SB 222.135E
SB 232.135E
SB 242.135E
SB 252.135E
SB 262.135E

SB 213.135E
SB 223.135E
SB 233.135E
SB 243.135E
SB 253.135E
SB 263.135E




150

outside
465 x 280 x 150 mm
inside
425 x 260 x 130 mm

SB 212.150E
SB 222.150E
SB 232.150E
SB 242.150E
SB 252.150E
SB 262.150E

SB 213.150E
SB 223.150E
SB 233.150E
SB 243.150E
SB 253.150E
SB 263.150E

Bottom non-perforated



Bottom perforated



3/4

New Safety Model

Variant with Quick Opening Mechanism

Please note that only the lid of the container is coloured.

SB 216.100E
SB 226.100E
SB 236.100E
SB 246.100E
SB 256.100E
SB 266.100E

SB 217.100E
SB 227.100E
SB 237.100E
SB 247.100E
SB 257.100E
SB 267.100E

100

outside
465 x 280 x 100 mm
inside
425 x 260 x 84 mm



SB 216.135E
SB 226.135E
SB 236.135E
SB 246.135E
SB 256.135E
SB 266.135E

SB 217.135E
SB 227.135E
SB 237.135E
SB 247.135E
SB 257.135E
SB 267.135E

135

outside
465 x 280 x 135 mm
inside
425 x 260 x 115 mm



SB 216.150E
SB 226.150E
SB 236.150E
SB 246.150E
SB 256.150E
SB 266.150E

SB 217.150E
SB 227.150E
SB 237.150E
SB 247.150E
SB 257.150E
SB 267.150E

150

outside
465 x 280 x 150 mm
inside
425 x 260 x 130 mm





Sterilization Container Systems

Bottoms and Lids

3/4

Bottom non-perforated



Article-No.	outside	inside
SB 205.100E	465 x 280 x 100 mm	425 x 262 x 84 mm
SB 205.135E	465 x 280 x 135 mm	425 x 262 x 121mm
SB 205.150E	465 x 280 x 150 mm	425 262 x 135 mm

silver
yellow
red
blue
green
black



465 x 280 mm

lid perforated



SB 200.010E
SB 200.020E
SB 200.030E
SB 200.040E
SB 200.050E
SB 200.060E

Safety Lid - Model - with exterior Press - Button
(Description on Page 5)

Bottom non-perforated



Article-No.	outside	inside
SB 206.100E	465 x 280 x 100 mm	425 x 262 x 85 mm
SB 206.135E	465 x 280 x 135 mm	425 x 262 x 121mm
SB 206.150E	465 x 280 x 150 mm	425 x 262 x 135 mm



465 x 280 mm

Safety Lid



Standard	New Model
SB 203.010	SB 203.010E
SB 203.020	SB 203.020E
SB 203.030	SB 203.030E
SB 203.040	SB 203.040E
SB 203.050	SB 203.050E
SB 203.060	SB 203.060E

3/4

Accessories

Article-No.	Dimensions
SB 290.030	405 x 255 x 30 mm
SB 290.050	405 x 255 x 50 mm
SB 290.070	405 x 255 x 70 mm
SB 290.100	405 x 255 x 100 mm




Stainless steel wire basket
intermediate crimp square hole wire mesh

Article-No.	Dimensions
SBS 290.030.02	405 x 255 x 35 mm
SBS 290.050.02	405 x 255 x 55 mm
SBS 290.070.02	405 x 255 x 75 mm
SBS 290.100.02	405 x 255 x 105 mm




Stainless steel screen
basket with lid

Article-No.	Dimensions
SBS 290.000	405 x 255 mm



Lid for stainless steel screen basket
with closure

Article-No.	Dimensions
SBS 290.030	405 x 255 x 30 mm
SBS 290.050	405 x 255 x 50 mm
SBS 290.070	405 x 255 x 70 mm
SBS 290.100	405 x 255 x 100 mm



Stainless steel screen basket




Sterilization Container Systems

Container Size

112

Standard Model



New Safety Model

Standard Model

Please note that only the lid of the container is coloured.


112

55

outside
285 x 280 x 55 mm
inside
270 x 270 x 33 mm



Bottom non-perforated



Bottom perforated



	SB 310.050E
	SB 320.050E
	SB 330.050E
	SB 340.050E
	SB 350.050E
	SB 360.050E

85

outside
285 x 280 x 85 mm
inside
270 x 270 x 64 mm



	SB 310.080E
	SB 320.080E
	SB 330.080E
	SB 340.080E
	SB 350.080E
	SB 360.080E

	SB 311.080E
	SB 321.080E
	SB 331.080E
	SB 341.080E
	SB 351.080E
	SB 361.080E

Bottom non-perforated



Article-No.	Outside	Inside
SB 305.050E	285 x 280 x 55 mm	270 x 270 x 33 mm
SB 305.080E	285 x 280 x 85 mm	270 x 270 x 64 mm



	SB 304.010E
	SB 304.020E
	SB 304.030E
	SB 304.040E
	SB 304.050E
	SB 304.060E

Lid perforated



Bottom perforated



Article-No.	Outside	Inside
SB 306.080E	285 x 280 x 85 mm	270 x 270 x 64 mm



	silver
	yellow
	red
	blue
	green
	black



Accessories

1/2


Article-No.	Dimensions
SB 390.030	255 x 240 x 55 mm
SB 390.050	255 x 240 x 85 mm

Stainless steel wire basket
intermediate crimp square hole
wire mesh



Article-No.	Dimensions
SBS 390.030.02	255 x 245 x 35 mm
SBS 390.050.02	255 x 245 x 55 mm

Stainless steel basket
with lid



Article-No.	Dimensions
SBS 390.000	255 x 245 mm

Lid for Stainless steel screen basket
with closure




Article-No.	Dimensions
SBS 390.030	255 x 245 x 30 mm
SBS 390.050	255 x 245 x 50 mm

Stainless steel screen basket



Bottom non-perforated



Bottom perforated



SB 310.100E
SB 320.100E
SB 330.100E
SB 340.100E
SB 350.100E
SB 360.100E

SB 311.100E
SB 321.100E
SB 331.100E
SB 341.100E
SB 351.100E
SB 361.100E

SB 310.135E
SB 320.135E
SB 330.135E
SB 340.135E
SB 350.135E
SB 360.135E

SB 311.135E
SB 321.135E
SB 331.135E
SB 341.135E
SB 351.135E
SB 361.135E

SB 310.150E
SB 320.150E
SB 330.150E
SB 340.150E
SB 350.150E
SB 360.150E

SB 311.150E
SB 321.150E
SB 331.150E
SB 341.150E
SB 351.150E
SB 361.150E

SB 310.200E
SB 320.200E
SB 330.200E
SB 340.200E
SB 350.200E
SB 360.200E

SB 311.200E
SB 321.200E
SB 331.200E
SB 341.200E
SB 351.200E
SB 361.200E

SB 310.260E
SB 320.260E
SB 330.260E
SB 340.260E
SB 350.260E
SB 360.260E

SB 311.260E
SB 321.260E
SB 331.260E
SB 341.260E
SB 351.260E
SB 360.260E

1/2

Standard Model

Please note that only the lid of the container is coloured.

100

outside
285 x 280 x 100 mm
inside
262 x 260 x 84 mm



135

outside
285 x 280 x 135 mm
inside
262 x 260 x 116 mm



150

outside
285 x 280 x 55 mm
inside
262 x 260 x 132 mm



200

outside
285 x 280 x 55 mm
inside
262 x 260 x 182 mm



260

outside
285 x 280 x 260 mm
inside
262 x 260 x 238 mm





Sterilization Container Systems

Standard Safety Model

Please note that only the lid of the container is coloured.


1/2




100

outside
285 x 280 x 100 mm
inside
262 x 260 x 84 mm

Bottom non-perforated



Bottom perforated



SB 312.100E
SB 322.100E
SB 332.100E
SB 342.100E
SB 352.100E
SB 362.100E

SB 313.100E
SB 323.100E
SB 333.100E
SB 343.100E
SB 353.100E
SB 363.100E



135

outside
285 x 280 x 135 mm
inside
262 x 260 x 116 mm

SB 312.135E
SB 322.135E
SB 332.135E
SB 342.135E
SB 352.135E
SB 362.135E

SB 313.135E
SB 323.135E
SB 333.135E
SB 343.135E
SB 353.135E
SB 363.135E



150

outside
285 x 280 x 150 mm
inside
262 x 260 x 132 mm

SB 312.150E
SB 322.150E
SB 332.150E
SB 342.150E
SB 352.150E
SB 362.150E

SB 313.150E
SB 323.150E
SB 333.150E
SB 343.150E
SB 353.150E
SB 363.150E



200

outside
285 x 280 x 200 mm
inside
262 x 260 x 182 mm

SB 312.200E
SB 322.200E
SB 332.200E
SB 342.200E
SB 352.200E
SB 362.200E

SB 313.200E
SB 323.200E
SB 333.200E
SB 343.200E
SB 353.200E
SB 363.200E




260

outside
285 x 280 x 260 mm
inside
262 x 260 x 238 mm


SB 312.260E
SB 322.260E
SB 332.260E
SB 342.260E
SB 352.260E
SB 362.260E

SB 313.260E
SB 323.260E
SB 333.260E
SB 343.260E
SB 353.260E
SB 363.260E

Bottom non-perforated



Bottom perforated



1/2

New Safety Model
Variant with Quick Opening Mechanism

Please note that only the lid of the container is coloured.

SB 316.100E
SB 326.100E
SB 336.100E
SB 346.100E
SB 356.100E
SB 366.100E

SB 317.100E
SB 327.100E
SB 337.100E
SB 347.100E
SB 357.100E
SB 367.100E

100

outside
285 x 280 x 100 mm
inside
262 x 260 x 84 mm



SB 316.135E
SB 326.135E
SB 336.135E
SB 346.135E
SB 356.135E
SB 366.135E

SB 317.135E
SB 327.135E
SB 337.135E
SB 347.135E
SB 357.135E
SB 367.135E

135

outside
285 x 280 x 135 mm
inside
262 x 260 x 116 mm



SB 316.150E
SB 326.150E
SB 336.150E
SB 346.150E
SB 356.150E
SB 366.150E

SB 317.150E
SB 327.150E
SB 337.150E
SB 347.150E
SB 357.150E
SB 367.150E

150

outside
285 x 280 x 150 mm
inside
262 x 260 x 132 mm



SB 316.200E
SB 326.200E
SB 336.200E
SB 346.200E
SB 356.200E
SB 366.200E

SB 317.200E
SB 327.200E
SB 337.200E
SB 347.200E
SB 357.200E
SB 367.200E

200

outside
285 x 280 x 200 mm
inside
262 x 260 x 182 mm



SB 316.260E
SB 326.260E
SB 336.260E
SB 346.260E
SB 356.260E
SB 366.260E

SB 317.260E
SB 327.260E
SB 337.260E
SB 347.260E
SB 357.260E
SB 367.260E

260

outside
285 x 280 x 260 mm
inside
262 x 260 x 238 mm





Sterilization Container Systems

Bottoms and Lids

1/2

Bottom non-perforated



Article-No.	outside	inside
SB 305.100E	285 x 280 x 100 mm	262 x 262 x 84 mm
SB 305.135E	285 x 280 x 135 mm	262 x 262 x 117 mm
SB 305.150E	285 x 280 x 150 mm	262 x 262 x 132 mm
SB 305.200E	285 x 280 x 200 mm	262 x 262 x 182 mm
SB 305.200E	285 x 280 x 260 mm	262 x 262 x 244 mm

Bottom perforated



Article-No.	outside	inside
SB 306.100E	285 x 280 x 100 mm	262 x 262 x 84 mm
SB 306.135E	285 x 280 x 135 mm	262 x 262 x 117 mm
SB 306.150E	285 x 280 x 150 mm	262 x 262 x 132 mm
SB 306.200E	285 x 280 x 200 mm	262 x 262 x 182 mm
SB 306.200E	285 x 280 x 260 mm	262 x 262 x 244 mm

silver
yellow
red
blue
green
black



lid perforated



SB 300.010E
SB 300.020E
SB 300.030E
SB 300.040E
SB 300.050E
SB 300.060E

Safety Lid- New Model- with exterior Press- Button (Description on Page 5)



285 x 280 mm



Safety Lid

Standard	New Model
SB 300.010E	SB 303.010E
SB 300.020E	SB 303.020E
SB 300.030E	SB 303.030E
SB 300.040E	SB 303.040E
SB 300.050E	SB 303.050E
SB 300.060E	SB 303.060E


1/2 Accessories

Article-No.	Dimensions
SB 390.030	255 x 245 x 30 mm
SB 390.050	255 x 245 x 50 mm
SB 390.070	255 x 245 x 70 mm
SB 390.100	255 x 245 x 100 mm



Stainless steel wire basket
intermediate crimp square hole wire mesh

Article-No.	Dimensions
SBS 390.030.02	255 x 245 x 35 mm
SBS 390.050.02	255 x 245 x 55 mm
SBS 390.070.02	255 x 245 x 70 mm
SBS 390.100.02	255 x 245 x 105 mm



Stainless steel screen basket
with lid

Article-No.	Dimensions
SBS 309.000	255 x 245 mm



Lid for Stainless steel screen basket
with closure

Article-No.	Dimensions
SBS 390.030	255 x 245 x 30 mm
SBS 390.050	255 x 245 x 50 mm
SBS 390.070	255 x 245 x 65 mm
SBS 390.100	255 x 245 x 100 mm



Stainless steel screen basket



Sterilization Container Systems

Mini Container

Please note that only the lid of the container is coloured.



40

outside
300 x 140 x 40 mm
inside
280 x 130 x 27 mm



70

outside
300 x 140 x 70 mm
inside
280 x 130 x 55 mm



100

outside
300 x 140 x 100 mm
inside
280 x 130 x 84 mm



Bottom non-perforated



Article-No.	outside	inside
SB 405.040E	300 x 140 x 40 mm	280 x 130 x 27 mm
SB 405.070E	300 x 140 x 70 mm	280 x 130 x 55 mm
SB 405.100E	300 x 140 x 100 mm	280 x 130 x 84 mm




Bottom perforated



Article-No.	outside	inside
SB 406.040E	300 x 140 x 40 mm	280 x 130 x 27 mm
SB R406.070E	300 x 140 x 70 mm	280 x 130 x 55 mm
SB R406.100E	300 x 140 x 100 mm	280 x 130 x 84 mm

Bottom non-perforated



SB 410.040E
SB 420.040E
SB 430.040E
SB 440.040E
SB 450.040E
SB 460.040E

Bottom perforated



SB 411.040E
SB 421.040E
SB 431.040E
SB 441.040E
SB 451.040E
SB 461.040E

SB 410.070E
SB 420.070E
SB 430.070E
SB 440.070E
SB 450.070E
SB 460.070E

SB 411.070E
SB 421.070E
SB 431.070E
SB 441.070E
SB 451.070E
SB 461.070E

SB 410.100E
SB 420.100E
SB 430.100E
SB 440.100E
SB 450.100E
SB 460.100E

SB 411.100E
SB 421.100E
SB 431.100E
SB 441.100E
SB 451.100E
SB 461.100E

silver
yellow
red
blue
green
black



Lid perforated



SB 400.010E
SB 400.020E
SB 400.030E
SB 400.040E
SB 400.050E
SB 400.060E

Mini Container Accessories

Article-No.	Dimensions
SB 490.050	265 x 120 x 45mm

Wire basket
for Mini container



Article-No.	Dimensions
SBS 490.050.02	271 x 117 x 45 mm

Stainless steel
screen basket for
Mini container
with lid



Article-No.	Dimensions
SBS 490.000	271 x 117 mm

Lid to stainless
steel screen basket
for Mini container



Article-No.	Dimensions
SBS 490.050	271 x 117 x 40 mm

Stainless steel
screen basket for
Mini container





Sterilization Container Systems

Dental Container Containers with Lids

Please note that only the lid of the container is coloured.

Bottom non-perforated



Bottom perforated



40

outside
310 x 190 x 40 mm
inside
280 x 182 x 35mm

SB 510.040E
SB 520.040E
SB 530.040E
SB 540.040E
SB 550.040E
SB 560.040E

SB 511.040E
SB 521.040E
SB 531.040E
SB 541.040E
SB 551.040E
SB 561.040E



65

outside
310 x 190 x 65 mm
inside
280 x 182 x 55 mm

SB 510.060E
SB 520.060E
SB 530.060E
SB 540.060E
SB 550.060E
SB 560.060E

SB 511.060E
SB 521.060E
SB 531.060E
SB 541.060E
SB 551.060E
SB 561.060E



80

outside
310 x 190 x 80 mm
inside
280 x 182 x 70 mm

SB 510.080E
SB 520.080E
SB 530.080E
SB 540.080E
SB 550.080E
SB 560.080E

SB 511.080E
SB 321.080E
SB 331.080E
SB 341.080E
SB 351.080E
SB 361.080E



100

outside
310 x 190x 100 mm
inside
280 x 182 x 90 mm

SB 510.100E
SB 520.100E
SB 530.100E
SB 540.100E
SB 550.100E
SB 560.100E

SB 511.100E
SB 521.100E
SB 531.100E
SB 541.100E
SB 551.100E
SB 561.100E



130

outside
310 x 190 x 130 mm
inside
280 x 182 x 124 mm

SB 510.130E
SB 520.130E
SB 530.130E
SB 540.130E
SB 550.130E
SB 560.130E

SB 511.130E
SB 521.130E
SB 531.130E
SB 541.130E
SB 551.130E
SB 561.130E

Dental Container Bottoms, Lids, and Accessories

silver
yellow
red
blue
green
black



Bottom perforated
without tag holders



Article-No.	outside	inside
SB 506.040E	310 x 190 x 40 mm	280 x 182 x 35 mm
SB 506.060E	310 x 190x 65 mm	280 x 182 x 55 mm
SB 506.080E	310 x 190 x 80 mm	280 x 182 x 70 mm
ER506.100E	310 x 190 x 100 mm	280 x 182 x 90 mm
SB 506.130E	310 x 190 x 130 mm	280 x 182 x 124 mm

Lid perforated
with tag holders



SB 500.010E
SB 500.020E
SB 500.030E
SB 500.040E
SB 500.050E
SB 500.060E



Bottom non-perforated
without tag holders



Article-No.	outside	inside
SB 505.040E	310 x 190 x 40 mm	280 x 182 x 35 mm
SB 505.060E	310 x 190x 65 mm	280 x 182 x 55 mm
SB 505.080E	310 x 190 x 80 mm	280 x 182 x 70 mm
SB 505.100E	310 x 190 x 100 mm	280 x 182 x 90 mm
SB 505.130E	310 x 190 x 130 mm	280 x 182 x 124 mm


Article-No.	Dimensions
SBS 590.050.02	265 x 175 x 43mm

Stainless steel screen basket also
with feet available

Article-No.	Dimensions
SBS 590.000	265 x 175 mm



Stainless steel screen
basket with lid



Lid to Stainless steel
screen basket



Dental Containers - Accessories



Wire basket
for Dental Cont

Article-No.	Dimensions
SB 590.045	265 x 170 x 45 mm



Screen basket

Article-No.	Dimensions
SBS 590.050	265 x 175 x 34 mm



Accessories

Article-No.	Dimensions
SB 800.000	190 0 mm

**Disposable Paper Filters
without Steam Indicator**
for 1/1, 3/4, 1/2, Container;
Pack of 100 pcs.



Article-No.	Dimensions
SB 802.000	235 x 118 mm

**Disposable Paper Filters
without Steam Indicator**
for Mini and Dental
Containers Pack of 100 pcs.



Article-No.	Color
SB 819.000	blue

Plastic Security Seal Pack of 100 pcs.



Accessories



Article-No. no marking	Article-No. with marking	Color
SB 820.010	SB 820.010B	gray
SB 820.020	SB 820.020B	yellow
SB 820.030	SB 820.030B	red
SB 820.040	SB 820.040B	blue
SB 820.050	SB 820.050B	green
SB 820.060	SB 820.060B	black
SB 820.070	SB 820.070B	violet
SB 820.080	SB 820.080B	orange
SB 820.100	SB 820.100B	brown

Aluminium Identifications Labels with Slot



Article-No. no marking	Article-No. with marking	Color
SB 821.020	SB 821.020B	gray
SB 821.030	SB 821.030B	yellow
SB 821.040	SB 821.040B	red
SB 821.050	SB 821.050B	blue
SB 821.060	SB 821.060B	green
SB 821.070	SB 821.070B	black
SB 821.080	SB 821.080B	violet
SB 821.100	SB 821.100B	brown

Aluminium Identifications Labels



Article-No. no marking	Article-No. with marking	Color
SB 822.010	SB 822.010B	gray
SB 822.020	SB 822.020B	yellow
SB 822.030	SB 822.030B	red
SB 822.040	SB 822.040B	blue
SB 822.050	SB 822.050B	green
SB 822.060	SB 822.060B	black
SB 822.070	SB 822.070B	brown

Aluminium Identifications Labels with Perforation



Strategic Business Holdings, Inc.

...Quality Products, Superior Customer Service, & Support